


Particulate Filter F1 HEPA H14

| References and Filter Type | |
|----------------------------|-----------|
| HP | H11011021 |

| Dimensions and Weight | |
|-----------------------|-------------------|
| Filter Dimensions | 190 x 403 x 61 mm |
| Protruding Seal(s) | 0 |
| Filter Weight | Env. 2.5 Kg |



| Filter Composition | |
|--------------------|----------------------|
| HEPA H14 Filter | Metal structure |
| | Non-woven fiber mesh |

| Features | |
|-----------------|---|
| HEPA H14 Filter | Particulate filtration for particles/powders or fumes |

| Filtration Efficiency | |
|------------------------|--|
| Particulate Filtration | 99.995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard |

| Warnings | |
|---|--|
| <p>Your protection is only guaranteed after approval by the Erlab® laboratory.</p> <p>Contact your supplier and benefit from the ESP® (Erlab Safety Program) to validate the compatibility of the filter type with your applications.</p> <ul style="list-style-type: none"> - This filter must be stored in its original packaging and in a dry place. - It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter. | |