


Particulate Filter GF4 HEPA H14

References and Filter Type	
HP	H11074031

Dimensions and Weight	
Filter Dimensions	745.5 x 394.5 x 89.5 mm
Protruding Seal(s)	0
Filter Weight	Env. 6 Kg



Filter Composition	
HEPA H14 Filter	Metal structure
	Non-woven fiber mesh

Features	
HEPA H14 Filter	Particulate filtration for particles/powders or fumes

Filtration Efficiency	
Particulate Filtration	99.995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard

Warnings	
Your protection is only guaranteed after approval by the Erlab® laboratory. Contact your supplier and benefit from the ESP® (Erlab Safety Program) to validate the compatibility of the filter type with your applications. - This filter must be stored in its original packaging and in a dry place. - It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.	