


## Particulate Filter HF4 ULPA U16

| References and Filter Type |           |
|----------------------------|-----------|
| U16                        | H11014061 |

| Dimensions and Weight |   |
|-----------------------|---|
| Filter Dimensions     | 757.5 x 406.5 x 89 mm   |
| Protruding Seal(s)    | 0   |
| Filter Weight         | <p><b>Warning:</b> Undefined array key "poids" in /var/www/clients/client0/web18/web/wp-content/themes/hello-elementor-child/pdf-templates/fiche_technique_filtre_hepa_ulpa_en.php on line 194</p> <p>Env. Kg</p> |



| Filter Composition |                      |
|--------------------|----------------------|
| ULPA U16 Filter    | Metal structure      |
|                    | Non-woven fiber mesh |

| Features        |   |
|-----------------|---|
| ULPA U16 Filter | Particulate filtration for particles/powders or fumes |

| Filtration Efficiency  |  |
|------------------------|--|
| Particulate Filtration | 99,99995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard |

## Warnings

**Your protection is only guaranteed after approval by the Erlab® laboratory.**

Contact your supplier and benefit from the **ESP® (Erlab Safety Program)** to validate the compatibility of the filter type with your applications.

- This filter must be stored in its original packaging and in a dry place.
- It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.



[www.erlab.com](http://www.erlab.com)

[info@erlab.net](mailto:info@erlab.net)