


Particulate Filter CF1 (lot de 2) HEPA H14

References and Filter Type	
HP	H12111021

Dimensions and Weight	
Filter Dimensions	190 x 403 x 61 mm
Protruding Seal(s)	0
Filter Weight	Env. 5 Kg



Filter Composition	
HEPA H14 Filter	Metal structure
	Non-woven fiber mesh

Features	
HEPA H14 Filter	Particulate filtration for particles/powders or fumes

Filtration Efficiency	
Particulate Filtration	99.995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard

Warnings	
<p>Your protection is only guaranteed after approval by the Erlab® laboratory.</p> <p>Contact your supplier and benefit from the ESP® (Erlab Safety Program) to validate the compatibility of the filter type with your applications.</p> <ul style="list-style-type: none"> - This filter must be stored in its original packaging and in a dry place. - It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter. 	