


Molecular Filters NR4301

References and Filter Types

AS	NR4301AS00	F	NR4301FF00
BE	NR4301BE00	K	NR4301KK00

Dimensions and Weight

Filter Dimensions	490 x 348 x 55 mm
Protruding Seal(s)	2
Filter Weight	Env. 9 Kg



Carbon Types

AS	Predominantly Organic Vapors
BE	Predominantly Acid Vapors
F	Predominantly Formaldehyde Vapors
K	Predominantly Ammonia Vapors

Filter Composition

Carbon Filter	Metal structure Granular activated carbon
---------------	--

Quality

Each filter is supplied with a quality control certificate indicating the filter's serial number, manufacturing date, and packaging date. The carbon used in the filter has been tested according to the following standards:

ASTM D2854-83	Method for determining the bulk density of activated carbon
ASTM D2862-82	Method for determining the particle size distribution of activated carbon
ASTM D5742-95	Method for determining the Butane working capacity of activated carbon
ASTM D2866-83	Method for determining the ash content of activated carbon
ASTM D2867-83	Method for determining the moisture content of activated carbon

Warnings

Your protection is only guaranteed after approval by the Erlab® laboratory.

Contact your supplier and benefit from the **ESP® (Erlab Safety Program)** to validate the compatibility of the filter type with your applications.

- This filter must be stored in its original packaging and in a dry place.
- It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.